

## **DONATED HUMAN TISSUE**

### **THIS ALLOGRAFT IS SUPPLIED STERILE**

This human tissue allograft is processed by CellRight® Technologies and distributed by Akorn Pharmaceuticals. All tissue was retrieved, processed, stored and distributed for use in accordance with the standards of the American Association of Tissue Banks (AATB), FDA requirements for Human Cellular and Tissue Based Products (HCT/Ps 21 CFR Part 1271), and applicable State regulations. The Donor has been determined to be eligible based on the results of screening and testing. Screening includes a review of medical and social history, available hospital records, infectious disease screening, autopsy report (if performed), and physical exam.

The Donor has been tested and was found negative or non-reactive for:

- Human Immunodeficiency Virus Types 1 and 2 Antibody (anti-HIV-1/anti-HIV-2)
- Hepatitis B Surface Antigen (HBsAg)
- Hepatitis B Core Antibody - Total (anti-HBc)
- Hepatitis C Virus Antibody (anti-HCV)
- Human Immunodeficiency Virus 1, Hepatitis C Virus, Hepatitis B Virus Nucleic Acid Test (HIV 1/HCV/HBV NAT)
- Syphilis Rapid Plasma Reagin or Treponemal Specific Assay
- West Nile Virus Nucleic Acid Test (WNV NAT)

Additional tests, including but not limited to HTLV I/II, may have been performed and were found to be acceptable for transplantation. U.S. Food and Drug Administration (FDA) licensed, approved, or cleared donor screening test kits are used when available. Communicable disease testing has been performed by a laboratory registered with the FDA to perform donor testing in accordance with the Clinical Laboratory Improvement Amendments of 1988 (CLIA) and 42 CFR part 493, or that has met equivalent requirements as determined by the Centers for Medicare and Medicaid Services (CMS). A list of additional communicable disease test(s) performed will be provided upon request.

CellRight Technologies Medical Director has determined this donor tissue to be suitable for transplantation. The testing and medical release records are maintained by CellRight Technologies. The names and addresses of the testing laboratories, the interpretation of all required infectious disease tests, and a listing of the documents reviewed as part of the relevant medical records are kept on file at CellRight Technologies and are available upon request.

Tissue has been processed using alcohol as a disinfectant. Tissue is air dried, not chemically dried. Tissue has been sterilized, using Cobalt 60, to a SAL of  $10^{-6}$  (Sterility Assurance Level).

### **WARNINGS AND PRECAUTIONS**

- Intended for use in one patient, on a single occasion only.
- Do not use if package integrity has been compromised. Once the user breaks the seal on the inner-most pouch, the tissue grafts must be transplanted or discarded.
- Tissue may not be sterilized or re-sterilized by your facility.
- This tissue is intended for use by qualified physicians.
- Although this tissue has been tested and screened for human pathogens, and processed under aseptic conditions, human derived tissue may still transmit infectious agents.
- It is the responsibility of the End-User clinician to maintain tissue intended for transplantation in appropriate storage conditions prior to transplant.

### **STORAGE**

AIR DRIED TISSUE - Maintain tissue at ambient temperature. Storage requirements are provided in this package insert and on the final label that accompanies each graft.

### **TISSUE PREPARATION**

#### **BEFORE USE – Examine Allograft Packaging – Do Not Use This Allograft If:**

1. Any of the package elements appear to be missing, tampered with or damaged.
2. The product label or identifying bar code is severely damaged, illegible or missing.
3. The expiration date shown on the package label has passed

If any of the above conditions exist or are suspected, this allograft should NOT be used.

**PREPARATION OF AIR-DRIED ALLOGRAFT TISSUE FOR USE**

1. Opening Peel Packages: peel outer package down and aseptically deliver inner peel pouch to the sterile field or sterile team member.
2. Remove tissue from Inner peel pouch.
  - a. Recommendation – cut the pouch close to the tissue, and gently separate the layers of the pouch in order to grasp the tissue with forceps.
3. Tissue should be applied dry; rehydration is not required.

**ADVERSE OUTCOMES**

Adverse outcomes potentially attributable to this tissue must be reported promptly to CellRight® Technologies.

Email: [CustomerCare@CellRightTechnologies.com](mailto:CustomerCare@CellRightTechnologies.com)

Phone: 1-210-659-9353

**TISSUE TRACKING**

Complete the enclosed Allograft Tracking Form and mail to CellRight Technologies. US Federal Regulations (21 CFR 1271.290(b)) requires proper tracking of this tissue. It is the responsibility of the end-user to provide this information, which enables CellRight Technologies to maintain records for the purpose of tracing the tissue post-transplant.

**Processed By:**

CellRight® Technologies  
1808 Universal City Blvd  
Universal City, TX 78148  
210-659-9353  
Fax: 210-659-9556

CellRight Technologies holds:

AATB Accreditation No. 00212  
US FDA Registration No. 3009234552  
Canadian Registration No. 100228  
California Tissue Bank ID No. CNC80949  
Florida License No. 212  
Maryland Tissue Bank No. TB1898  
New York State Tissue Bank ID No. 1779

CellRight Technologies is a Registered Tissue Bank in the following state(s):

Oregon  
Delaware

**Distributed By:**



Akorn Pharmaceuticals  
1925 West Field Court, Suite 300  
Lake Forest, IL 60045  
[www.Akorn.com](http://www.Akorn.com)  
[www.AcellFX.com](http://www.AcellFX.com)